



General Specifications				
LABORATORY FREEZERS TO -30°C				
MODEL NO.	CU. FT /LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS, INCHES, (CM) W x F-B x H	EXTERIOR DIMENSIONS, INCHES (CM) W x F-B x H
<a href="#">DHF29-23SD</a>	23.1 / 654	-15°C to -30°C	22" x 28" x 60" (56 x 71 x 152)	26" x 33.5" x 84.5" (66 x 85 x 215)
<a href="#">DHF29-46SD</a>	46.6 / 1320	-15°C to -30°C	48" x 28" x 60" (122 x 71 x 152)	52" x 33.5" x 84.5" (132 x 85 x 215)
<a href="#">DHF29-74SD</a>	74 / 2095	-15°C to -30°C	74" x 28" x 60" (188 x 71 x 152)	78" x 33.5" x 84.5" (198 x 85 x 215)
<a href="#">DHF30-30MDS</a>	30 / 849	-10°C to -25°C	32" x 22.5" x 65" (81 x 57 x 165)	35.5" x 35" x 80" (90 x 89 x 203)
<a href="#">DHN20-24MDS</a>	24	-10°C to -25°C	23.5" x 30" x 59" (60 x 74 x 150)	27.5" x 34.88" x 81.5" (70 x 89 x 207)
<a href="#">DHW20-20MDP</a>	19.7 / 558	-20°C	28" x 16.25" x 56"	33" x 27.5" x 76.5"
REFRIGERATOR / FREEZER COMBINATION UNITS				
<a href="#">DHF-GPRF</a>	17.9 / 509	2°C to 5°C (R) -18°C (F)		30" x 30.75" x 66"
<a href="#">DHC-20RF</a>	20 / 560	4°C (R) / -20°C (F)		26" x 35.5" x 82.25"
<a href="#">DHC-40RF</a>	40	4°C (R) / -20°C (F)		52" x 35.5" x 83.25"
GLASS DOOR LABORATORY FREEZERS TO -22°C				

<a href="#">DHC20-20GD</a>	20 / 566	-12°C to -22°C		26"x 35.5"x 86.75" (66 x 90 x 220)
<a href="#">DHC20-48GD</a>	48 / 1360	-12°C to -22°C		52"x 35.5"x 85.75" (132 x 90 x 218)
<a href="#">DHC20-70GD</a>	70 / 1980	-12°C to -22°C		78"x 35.5"x 85.75" (198 x 90 x 218)