

# Laboratory furnaces L for asphalt test

Electric furnace SNOL 45/700L is designed for test with asphalt, bitumen, heat treatment of metal, ceramic up to 700°C.

Heating elements on the tube

Three sided heating

- Chamber build from refractory fibre and bricks
- A tray with a lid and grille
- Two gaps for fresh air inlet are in the underpart of door
- The tube of air outlet is on the top of furnace. Tube is connected with air ventilator
- Microprocessor temperature controller
- High quality thermoinsulation materials
- Door interlock switch
- Ceramic bottom plate
- High precision level
- Available connection to PC
- Two years limited warranty and later service

Type	T max C	Chamber dimensions, mm			Volume L	Outside dimensions, mm			Power kW	Voltage V	Weight kg
		W	D	H		W	D	H			
SNOL 45/700 L	700	400	450	250	45	810	955	1610	9	400	280