

# Shaft furnaces - LSh

Universal precise laboratory shaft furnace is suitable for many applications in the laboratory.

- T max. 900°C
- Heating is from four sides
- Heating elements are embedded in the grooves
- Micro processor temperature controller
- High quality thermoinsulation materials
- Ceramic chamber
- Door interlock switch
- Ceramic bottom plate
- Hermetic closing door
- High precision level
- Available connection to PC
- Two years limited warranty and latter service

Type	T max C	Chamber dimensions, mm			Volume L	Outside dimensions, mm			Power kW	Voltage V	Weight kg
		W	D	H		W	D	H			
<u>SNOL</u> <u>10/900</u> <u>LSh</u>	900	150	150	450	<b>10</b>	860	750	800	4,5	230	100