

# Tube furnace - LV

Standard laboratory electric tube furnace can be operated up to 1250°C. This model could be used universally for many processes.

- Ceramic tube working chamber
- Micro processor temperature controller
- High quality thermoinsulation materials
- Low power consumption
- Two years limited warranty and latter service
- Available connection to PC

Type	T max C	Chamber dimensions, mm			Volume L	Outside dimensions, mm			Power kW	Voltage V	Weight kg
		W	D	H		W	D	H			
<u>SNOL</u> 0,2/1250 <u>LV</u>	1250	Ø 25		100	0,2	470	340	480	1,7	230	19