



General Specifications TO -40° C (-40°F)				
CHEST, HORIZONTAL FREEZERS				
MODEL NO.	CU. FT /LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS, INCHES, (CM) W x F-B x H	EXTERIOR DIMENSIONS, INCHES (CM) W x F-B x H
CH45-5	5 - 142	0° C to -45° C	28 1/2" x 16 1/2" x 22" (72.4 x 41.9 x 55.9)	38 1/2" x 24 1/2" x 27 1/2" (97.8 x 62.2 x 69.9)
CH43-9	9 - 255	0° C to -43° C	28 1/2" x 16 1/2" x 36" (72.4 x 41.9 x 91.4)	38 1/2" x 24 1/2" x 42" (97.8 x 62.2 x 106.7)
CH40-13	13 - 370	0° C to -40° C	28 1/2" x 22" x 37 1/2" (72.4 x 55.9 x 95.3)	38 1/2" x 31" x 43" (97.8 x 78.7 x 109.2)
UPRIGHT, VERTICAL FREEZERS				
U40-13	13 - 370	0° C to -40° C	20" x 22" x 51" (51 x 56 x 130)	35" x 33" x 79" (89 x 84 x 201)
U40-18	18 - 505	0° C to -40° C	23" x 26" x 51" (58 x 66 x 130)	36" x 35 1/2" x 79" (90 x 90 x 201)
U40-22	22 - 626	0° C to -40° C	28 1/2" x 26" x 51" (72 x 66 x 130)	41 1/2" x 35 1/2" x 79" (105 x 90 x 201)
U40-25	25 - 710	0° C to -40° C	34" x 25" x 51- 1/4" (86.4 x 63.5 x 130)	49" x 35-1/2" x 79" (124.5 x 90.2 x 201)

